**Seago SG-972-S5**

-

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## **Představení**

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**User's guide**

## **Performance**

The Seago S5 is a sonic toothbrush with five different cleaning modes, 40,000 oscillations per minute, IPX7 durability and USB recharging.

## **Contents of the package**

Toothbrush, travel case, two heads for brushing teeth, charging cable, instructions.

## **First use**

Charge the toothbrush. Attach the brush head and make sure it is firmly connected to the base. Apply a reasonable amount of paste to the toothbrush. Attach the toothbrush to the teeth and turn it on. The toothbrush is adjusted to always give a brief interruption of vibration after 30 seconds to the point that they should change the cleaning area. The total cleaning time is set to 2 minutes. Then the toothbrush turns off.

## **Cleaning modes**

Use the multifunction button to switch cleaning modes. By repeatedly pressing the button, you switch between individual cleaning modes (sorted from above) – bleaching, cleaning, sensitive

cleaning, polishing, massage. The selected mode is indicated by the diode on.

## **Recharge**

Low battery capacity indicates a flashing red light. Connect the cable with the connector to the hole at the bottom of the brush. Open the cap and connect the cable. When charging, the LED on the brush turns red. If the battery is low for a long time, it may take longer than led charging is indicated.

## **Deployment / replacement of the head**

Brush heads can only be fitted to the body in one way where the brush bristles are parallel to the multifunction button. You can also see the correct deployment on the sales packaging of the product. Gently place the head on the metal protrusion of the body. If the head is in the right direction, it will push you all the way to the brush body. If it's deployed obliquely, it won't push to the body.

## **Safety instructions**

The toothbrush is intended only for brushing teeth.

If you have recently had surgery in your oral cavity, talk to your dentist about using your toothbrush.

If your gums don't stop bleeding while you're in use, talk to your dentist about using your toothbrush.

If children use the toothbrush, they should be checked at least at the beginning.

## Disposal

The product contains batteries and/or recyclable electrical waste. In order to protect the environment, do not dispose of the product with ordinary waste, but

hand it over for recycling to a collection point of electrical waste.

## Declaration of conformity

The importer, Satomar, s.r.o. (Karlova 37, 614 00, Brno) declares that the product meets all the essential requirements of eu

directives that apply to it. The full text of the Declaration of Conformity can be downloaded <http://www.satomar.cz> . Copyright © 2021

Satomar, s.r.o. All rights reserved.

Appearance and specifications may be changed without notice. Print errors reserved. The most up-to-date version of the manual can always be found [www.mobilnipodpora.cz](http://www.mobilnipodpora.cz)

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